


Zimmer Biomet Instrument Kit System

Reference Guide

Tapered Screw-Vent® (TSV™) Implant System
Trabecular Metal™ Dental Implants
3.1 mmD Eztetic® Dental Implants



Staging Block




NP Surgical Module for
3.1 mmD Eztetic Implants



ZIMMER BIOMET
Your progress. Our promise.®

Instrument Kit System For TSV, Trabecular Metal And 3.1 mmD Eztetic Implants






























TSV Surgical Kit
(TSVKIT)

- | | | | | | | | | | | |
|--|---|--------------------------------------|--|---|---|---|--|---|---|---|
| 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | 10 | 11 |
| 3.0 mmD
Round Bur
1203 | 2.3 mmD
Drill,
22 mmL
SV2.3DN | | 2.8 mmD
Drill,
22 mmL
SV2.8DN | 3.4
2.8 mmD
Step Drill,
22 mmL
TSV3DN | 3.4 mmD
Drill,
22 mmL
SV3.4DN | 3.8/
3.4 mmD
Step Drill,
22 mmL
TSV3.8DN | 3.8 mmD
Drill,
22 mmL
SV3.8DN | 4.4/
3.8 mmD
Step Drill,
22 mmL
TSV4DN | 5.1 mmD
Drill,
22 mmL
SV5.1DN | 5.7/5.1
mmD
Step Drill,
22 mmL
TSV6DN |
| 12 | 13 | 14 | 15 | 16 | 17 | 18 | 19 | 20 | 21 | |
| Tapered
Pilot Drill,
2.1/
1.6 mmD,
8 mmL
0201 | 2.3 mmD
Drill,
16 mmL
SV2.3DSN | Paralleling
Tool (Qty: 4)
PPAR | 2.8 mmD
Drill,
22 mmL
SV2.8DN | 3.4/
2.8 mmD
Step Drill,
16 mmL
TSV3DSN | 3.4 mmD
Drill,
16 mmL
SV3.4DSN | 3.8/
3.4 mmD
Step Drill,
16 mmL
TSV3.8DSN | 3.8 mmD
Drill,
16 mmL
SV3.8DSN | 4.4/
3.8 mmD
Step Drill,
16 mmL
TSV4DSN | 5.1 mmD
Drill,
16 mmL
SV5.1DSN | 5.7/
5.1 mmD
Step Drill,
16 mmL
TSV6DSN |
| 22 | 23 | 24 | 25 | 26 | 27 | 28 | 29 | | | |
| Drill
Extender
DE | 3.7 mmD
Bone Tap
TT3.7 | 4.1 mmD
Bone Tap
TT4.1 | 4.7 mmD
Bone Tap
TT4.7 | 6.0 mmD
Bone Tap
TT6.0 | 2.5 mm GemLock
Hex Drill
RHD2.5 | 2.5 mm GemLock
Hex Tool, Short
RH2.5 | 2.5 mm GemLock
Hex Tool, Long
RHL2.5 | | | |

The NP Surgical Module snaps into the Tapered Screw-Vent Surgical Kit



NP Surgical Module (NPMOD)

							
2.4 mmD Driva™ Drill, 22 mmL SV2.4DN	2.8/2.4 mmD Driva Surgical Step Drill, 22 mmL ZOP28DN	3.0 mmD Cortical Bone Tap Tool ZOPTT30	2.1 mm GemLock® Hex Drill, Latch Lock, Short CHD2.1	2.1 mm GemLock Hex Drill, Latch Lock, Long CHDL2.1	2.1 mm GemLock Hex Tool, Ratchet, Short CHR2.1	2.1 mm GemLock Hex Tool, Ratchet, Long CHRL2.1	Removal Tool for Eztetic Abutments CLRT2
							
Drill Extender DE	3.7 mmD Bone Tap TT3.7	4.1 mmD Bone Tap TT4.1	4.7 mmD Bone Tap TT4.7	6.0 mmD Bone Tap TT6.0	2.5 mm GemLock Hex Drill RHD2.5	2.5 mm GemLock Hex Tool, Short RH2.5	2.5 mm GemLock Hex Tool, Long RHL2.5
							
3.0 mm Hex Insertion Drill HX3.0D	3.0 mm Hex Insertion Tool, 17 mmL HX3.0-S	3.0 mm Hex Insertion Tool, 25 mmL HXL3.0-S	Removal Tool TLRT2	1.25 mm Hex Driver with GemLock Retention, 22 mmL HXGR1.25	1.25 mm Hex Driver with GemLock Retention, 30 mmL HXLGR1.25	1.25 mm Hex Tool, 17 mmL HX1.25	1.25 mm Long Hex Tool, 22 mmL HXL1.25
							
1.25 mm Hex Drill 23 mmL HX1.25D	GemLock Retaining Square Ratchet RSR	Screwdriver Handle with Square Connection SSHS					

Cleaning of Instruments*

1. Disassemble multi-piece components.
2. Rinse instruments in cool to lukewarm drinkable water for 2 and a 1/2 minutes.
3. For drills, use the cleaning wire to remove any debris from the irrigation channel. Using a 25 gauge needle, flush the drill lumen with water to remove any remaining debris.
4. Sonicate the instruments for 10 minutes in an ultrasonic cleaner with a pH-neutral enzymatic detergent diluted with tap water per the manufacturer's instructions.
5. Rinse the instruments with drinkable tap water for 3 minutes.
6. Inspect the instruments for signs of wear, damage, or unrecognizable color identification and replace the instruments accordingly.

Cleaning of Trays and Staging Block*

1. Remove all parts and insert from the surgical tray. Clean parts per above instructions.
2. Rinse the tray and tray insert with cool to lukewarm drinkable tap water to remove all visible soil.
3. Fully immerse the kit in enzymatic detergent, prepared per manufacturer's specifications, and allow the kit to soak for a minimum of one minute.
4. Use a damp cloth or soft-bristle brush to wipe and remove any excess soil from each part.
5. Rinse thoroughly with tap water for 3 minutes.
6. Dry components. Reassemble kit contents and follow sterilization guidelines.

*For detailed cleaning and sterilization instructions, refer to the Instructions for Use provided with the products.

Sterilization*

1. Individual parts should be placed in a sterilization pouch prior to sterilization.
2. Kits should be populated with clean instruments, placed in a sterilization pouch and sealed.
3. Validated sterilization parameters:

Cycle Type	Temperature	Exposure Time	Dry Time
Gravity (steam)	132°C 270°F	15 mins	20 mins
Pre-vacuum (steam)	132°C 270°F	4 mins	20 mins
Pre-vacuum (steam)	134°C 273°F	3 mins	20 mins
Pre-vacuum (steam)	134°C 273°F	18 mins	20 mins




For maximum cutting efficiency, replace drills frequently.

Tapered Screw-Vent And Trabecular Metal Implant Surgical Protocol

Intuitive Flow And Color-Coding


A simple color-coding system identifies drills for each implant diameter, allowing you to easily follow any surgical sequence step-by-step. As an example, surgical drills required for placement of the 3.7 mmD Tapered Screw-Vent Implant are represented by horizontal green bars on the kit surface and are logically organized in the order you would use them from left to right. The color-coding also allows you to easily identify your drill options for soft- or dense-bone protocols – a dotted color bar denotes a final soft-bone drill, while the following solid color bar denotes a final dense-bone drill.

Tapered Screw-Vent And Trabecular Metal Implants Color-Coding



Step 1

The 3.7 mmD Tapered Screw-Vent and Trabecular Metal Implant are color-coded in green. Start with the first green bar on the kit, which indicates the first drill to be used in the drilling sequence for this implant size.



Step 2

Follow the green color bars from left to right. In a soft-bone protocol, the dotted green bar represents the final drill. For dense bone, skip the dotted green bar and move on directly to the next solid green bar. The last solid bar in the sequence represents the final drill for dense bone.




Step 3

When drilling in dense bone, you can optionally use the 3.7 mmD cortical bone tap located in a green grommet directly below the last solid green bar in the sequence.


3.7 mmD Tapered Screw-Vent and Trabecular Metal Implant (3.5 mmD Platform)




3.7 mmD




1
SV2.3DN
2.3 mmD
Drill



2
FOR SOFT BONE
SV2.8DN
2.8 mmD
Drill




2
FOR DENSE BONE
TSV3DN
3.4/2.8 mmD
Drill



3
OPTIONAL FOR
DENSE BONE
TT3.7
3.7 mmD
Cortical Bone Tap




4.1 mmD Tapered Screw-Vent and Trabecular Metal Implant (3.5 mmD Platform)




*When placing the 4.1 mmD Trabecular Metal Dental Implant in dense bone (Type D1), add an additional drill step utilizing the SV3.8DN/SV3.8DSN drill after TSV3.8DN/TSV3.8DSN.

4.7 mmD Tapered Screw-Vent and Trabecular Metal Implant (4.5 mmD Platform)



6.0 mmD Tapered Screw-Vent and Trabecular Metal Implant (5.7 mmD Platform)



*In dense bone, an optional additional step drill may be used before TSV6DN/TSV6DSN: TSV5.1DN/TSV5.1DSN. Note this additional drill is sold separately and is not included in kits.

Eztetic Implant Surgical Protocol

Intuitive Flow And Color-Coding


A simple color-coding system identifies drills for each implant diameter, allowing you to easily follow any surgical sequence step-by-step. As an example, surgical drills required for placement of the 3.1 mmD Eztetic Implant are represented by horizontal blue bars on the kit and NP Module and are logically organized in the order you would use them from left to right and up to the NP Module. The color-coding also allows you to easily identify your drill options for soft- or dense-bone protocols – a dotted color bar denotes a final soft-bone drill, while the following solid color bar denotes a final dense-bone drill.



3.1 mmD Eztetic Implant Color-Coding


With the instrument kit system comes a simple way of working. Its unique, color-coded surgical protocol labeling system helps to guide you effortlessly through each drilling sequence.

3.1 mm Implant Diameter



Step 1

The 3.1 mmD Eztetic Implant is color-coded in blue. Start with the first blue bar in the main kit.



Step 2


Follow the blue color bars from left to right for the 2.3 mmD Drill, and up to the NP Surgical Module. Utilize the 2.4 mmD Drill as the final drill in soft bone. The last blue solid bar in the sequence is located in the NP Surgical Module and represents the final drill in dense bone (2.8/2.4 mmD).




Step 3

When drilling in dense bone, you can optionally use the 3.0 mmD cortical bone tap located in the blue grommet that follows the 2.8/2.4 mm drill.


3.1 mmD Eztetic Implants




0201DSN
2.1/1.6 mmD,
8 mmL - 11.5 mmL
Drill




SV2.3DN
2.3 mmD
Drill




SV2.4DN
2.4 mmD
Drill



ZOP28DN
2.8/2.4 mmD
Drill



**OPTIONAL FOR
DENSE BONE
ZOPTT30**
3.0 mmD
Cortical Bone Tap





Contact us at 1-800-342-5454 or visit
zimmerbiometdental.com

Zimmer Biomet Dental
Global Headquarters
4555 Riverside Drive
Palm Beach Gardens, FL 33410
Phone: +1-561-776-6700
Fax: +1-561-776-1272

Unless otherwise indicated, as referenced herein, all trademarks are the property of Zimmer Biomet; and all products are manufactured by one or more of the dental subsidiaries of Zimmer Biomet Holdings, Inc. and marketed and distributed by Zimmer Biomet Dental and its authorized marketing partners. For additional product information, please refer to the individual product labeling or instructions for use. Product clearance and availability may be limited to certain countries/regions. This material is intended for clinicians only and does not comprise medical advice or recommendations. Distribution to any other recipient is prohibited. This material may not be copied or reprinted without the express written consent of Zimmer Biomet Dental. ZB0105 REV D 03/21 ©2021 Zimmer Biomet. All rights reserved.

