

# TSV™ BellaTek® Laboratory Solutions

Optimization By Design®



# Digital Dentistry Solutions for TSV Implants

## Business Optimization With TSV™ BellaTek® Digital Dentistry Solutions

Optimize the solutions you give to your customers with BellaTek digital dentistry solutions. Zimmer Biomet Dental offers your laboratory a full range of options for aesthetic restorations with a knowledgeable support team to help you reach your business objectives.



## Aesthetic Outcomes

Fully digital workflows and laboratory designed products are patient specific, allowing you more control over restorations for better aesthetic and functional outcomes compared to traditional procedures.



## Customized Treatment Planning Options

Let Zimmer Biomet Dental help you focus on providing your customers with highly aesthetic outcomes for their patients with patient-specific BellaTek Abutments for TSV Implants.

## Efficiency

The Zimmer Biomet Dental Team helps to make your business more efficient with exceptional support from sales, customer service and technical support teams with best-in-class service and assistance to improve restorative options for clinicians.



## Customized Product Solutions Offering Aesthetic Patient Outcomes

### BellaTek Impression System



- A unique and proprietary, simplified process, which is designed to open up new workflows to laboratories. The result is a highly aesthetic BellaTek Definitive Patient-Specific Abutment for TSV Implants.
- There is no need to create a cast, which results in less steps in the treatment process compared to traditional procedures.

### TSV Implant Digital Laboratory Analogs



- The TSV Implant digital lab analogs are designed to be used with 3D printed models to support the fabrication of definitive restorations in the digital workflow. The analogs replicate the position and orientation of the TSV Implant in the printed model.
- The TSV Implant digital lab analogs incorporate a custom-designed feature to protect the TSV friction-fit connection during fabrication and prior to placement of the final restoration. A Lab Analog Placement Tool enables the analog to be installed into the model with ease and accuracy. One-position locking design: proprietary anti-rotation feature to ensure accurate hex timing within the model.



### Encode Empowered Laboratory (EEL) Designed BellaTek Abutments

All design work is completed by trusted certified laboratory technicians at your local EEL.

The Laboratory Designed Abutment is a solution providing dental laboratories with more control over abutment design coupled with the quality expected from Zimmer Biomet Dental. Additionally, there are reduced shipping costs and quicker turnaround times compared to traditional procedures.





## Become a BellaTek® Encode® Empowered Laboratory

Through the development of the BellaTek Encode Empowered program, the BellaTek Encode Impression System becomes even more powerful in the hands of laboratories around the world!

Due to advancements in digital dentistry, BellaTek Laboratories are now able to control the design of their CAD/CAM BellaTek Encode Abutments, experience faster turnaround times as well as enjoy a more simplified workflow to better support and serve their clinicians. Through the development of this program, laboratories are provided with the following benefits:

- Scan and design control from a BellaTek Encode Impression
- Reduced turnaround times
- Abutment material flexibility
- Cost effective milling options

This is accomplished through a non-exclusive license of the BellaTek Encode Library to participating laboratories throughout the USA and Canada, resulting in a timelier and more cost efficient workflow for both the laboratory and restoring clinician.

Please contact your local representative for more information on becoming a BellaTek Encode Empowered Laboratory or speak directly to our dedicated BellaTek Customer Service Center at 1-888-800-8045.

*Bella***Tek**®  
*Encode*® *Empowered Laboratory*



Contact us at 1-800-342-5454 or visit  
[zimmerbiometdental.com](http://zimmerbiometdental.com)

Zimmer Biomet Dental  
Global Headquarters  
4555 Riverside Drive  
Palm Beach Gardens, FL 33410  
Tel: +1-561-776-6700  
Fax: +1-561-776-1272

Unless otherwise indicated, as referenced herein, all trademarks are the property of Zimmer Biomet; and all products are manufactured by one or more of the dental subsidiaries of Zimmer Biomet Holdings, Inc. and marketed and distributed by Zimmer Biomet Dental and its authorized marketing partners. For additional product information, please refer to the individual product labeling or instructions for use. Product clearance and availability may be limited to certain countries/regions. This material is intended for clinicians only and does not comprise medical advice or recommendations. Distribution to any other recipient is prohibited. This material may not be copied or reprinted without the express written consent of Zimmer Biomet Dental.  
ZB0522 REV B 12/19 ©2019 Zimmer Biomet. All rights reserved.

