

Screw Removal Tools Kit

Certain® And External Hex Implant Systems



The Right Tools For A Rare Occasion

1 Remove And Retrieve

A restorative solution designed to extract broken abutment screws in the unlikely event of screw fracture.

2 Design Precision

Increase the odds of rescuing the fractured screw with tools specifically designed and manufactured for Certain and External Hex Implant Systems.

3 Convenience

Each kit contains all necessary parts to facilitate screw removal, or tools can be purchased individually.



ZIMMER BIOMET
Your progress. Our promise.®

Drill Guide Handle

Aligns the reversal drill over the broken screw inside of the implant



Screw Removal Extraction Tool

Captures and removes the broken screw



Waxing Screw/Guide Pin

Verifies the internal threads of the implant



Manual Reversing Drill

Reverses the broken screw from the internal threads of both Certain and External Hex implants.



Ordering Information

ISRT10N - Certain Components Kit

ISRT05N	ISRT06N	USRT07	IWSU30
Certain Drill Guide Handle	Certain Screw Removal Extraction Tool	Manual Reversing Drill	Certain Waxing Screw/Guide Pin

SRT10N - External Hex Components Kit

SRT05N	SRT06N	USRT07	WSU30
External Hex Drill Guide Handle	External Hex Screw Removal Extraction Tool	Manual Reversing Drill	External Hex Waxing Screw/Guide Pin

Please refer to ZBINST0021 for Recommended Guidelines and Instructions for Use.

Contact us at 1-800-342-5454 or visit zimmerbiometdental.com

Zimmer Biomet Dental
Global Headquarters
4555 Riverside Drive
Palm Beach Gardens, FL 33410
Phone: +1-561-776-6700
Fax: +1-561-776-1272

Unless otherwise indicated, as referenced herein, all trademarks are the property of Zimmer Biomet; and all products are manufactured by one or more of the dental subsidiaries of Zimmer Biomet Holdings, Inc. and marketed and distributed by Zimmer Biomet Dental and its authorized marketing partners. For additional product information, please refer to the individual product labeling or instructions for use. Product clearance and availability may be limited to certain countries/regions. This material is intended for clinicians only and does not comprise medical advice or recommendations. Distribution to any other recipient is prohibited. This material may not be copied or reprinted without the express written consent of Zimmer Biomet Dental. ZB0617 REV C 12/19 ©2019 Zimmer Biomet. All rights reserved.

