

Prosthetic Overview









Trabecular Metal™, TSV® and Eztetic®
Implant Systems

















ZIMMER BIOMET


































Your progress. Our promise.®

Trabecular Metal, Tapered Screw-Vent® and Eztetic Implants

Implants		Trabecular Metal			
		Ø 3.7 mm diameter	Ø 4.1 mm diameter	Ø 4.7 mm diameter	Ø 6.0 mm diameter
Machined Collar (0.5 mm) With Microgrooves		 TMMB10, 10 mm TMMB11, 11.5 mm TMMB13, 13 mm TMMB16, 16 mm	 TMM4B10, 10 mm TMM4B11, 11.5 mm TMM4B13, 13 mm	 TMMWB10, 10 mm TMMWB11, 11.5 mm TMMWB13, 13 mm	 TMM6B10, 10 mm TMM6B11, 11.5 mm TMM6B13, 13 mm
	MTX Textured Collar With Microgrooves	 TMTB10, 10 mm TMTB11, 11.5 mm TMTB13, 13 mm TMTB16, 16 mm	 TMT4B10, 10 mm TMT4B11, 11.5 mm TMT4B13, 13 mm	 TMTWB10, 10 mm TMTWB11, 11.5 mm TMTWB13, 13 mm	 TMT6B10, 10 mm TMT6B11, 11.5 mm TMT6B13, 13 mm














Implants		Eztetic	Tapered Screw-Vent			
		Ø 3.1 mm diameter	Ø 3.7 mm diameter	Ø 4.1 mm diameter	Ø 4.7 mm diameter	Ø 6.0 mm diameter
Machined Collar (1 mm) Without Microgrooves			 TSVB8, 8 mm TSVB10, 10 mm TSVB11, 11.5 mm TSVB13, 13 mm TSVB16, 16 mm	 TSV4B8, 8 mm TSV4B10, 10 mm TSV4B11, 11.5 mm TSV4B13, 13 mm TSV4B16, 16 mm	 TSVWB8, 8 mm TSVWB10, 10 mm TSVWB11, 11.5 mm TSVWB13, 13 mm TSVWB16, 16 mm	 TSV6B8, 8 mm TSV6B10, 10 mm TSV6B11, 11.5 mm TSV6B13, 13 mm TSV6B16, 16 mm
	Machined Collar (0.5 mm) With Microgrooves	 CM318, 8 mm CM3110, 10 mm CM3111, 11.5 mm CM3113, 13 mm CM3116, 16 mm	 TSVMB8, 8 mm TSVMB10, 10 mm TSVMB11, 11.5 mm TSVMB13, 13 mm TSVMB16, 16 mm	 TSM4B8, 8 mm TSM4B10, 10 mm TSM4B11, 11.5 mm TSM4B13, 13 mm TSM4B16, 16 mm	 TSMWB8, 8 mm TSMWB10, 10 mm TSMWB11, 11.5 mm TSMWB13, 13 mm TSMWB16, 16 mm	 TSM6B8, 8 mm TSM6B10, 10 mm TSM6B11, 11.5 mm TSM6B13, 13 mm TSM6B16, 16 mm
MTX Textured Collar With Microgrooves	 CT318, 8 mm CT3110, 10 mm CT3111, 11.5 mm CT3113, 13 mm CT3116, 16 mm	 TSVTB8, 8 mm TSVTB10, 10 mm TSVTB11, 11.5 mm TSVTB13, 13 mm TSVTB16, 16 mm	 TSVT4B8, 8 mm TSVT4B10, 10 mm TSVT4B11, 11.5 mm TSVT4B13, 13 mm TSVT4B16, 16 mm	 TSVTWB8, 8 mm TSVTWB10, 10 mm TSVTWB11, 11.5 mm TSVTWB13, 13 mm TSVTWB16, 16 mm	 TSVT6B8, 8 mm TSVT6B10, 10 mm TSVT6B11, 11.5 mm TSVT6B13, 13 mm TSVT6B16, 16 mm	

Prosthetic Components for Trabecular Metal, Tapered Screw-Vent and Eztetic Implants

Healing Abutments	Eztetic Implants		Trabecular Metal and Tapered Screw-Vent Implants									
	Ø 2.9 mm platform		Ø 3.5 mm platform		Ø 4.5 mm platform		Ø 5.7 mm platform					
Encode Healing Abutments	N/A		 TEHA3383, Ø 3.8 mm/3 mm TEHA3385, Ø 3.8 mm/5 mm TEHA3387, Ø 3.8 mm/7 mm TEHA3503, Ø 5 mm/3 mm TEHA3505, Ø 5 mm/5 mm	 TEHA4503, Ø 5.0 mm/3 mm TEHA4505, Ø 5.0 mm/5 mm TEHA4563, Ø 5.6 mm/3 mm TEHA4565, Ø 5.6 mm/5 mm TEHA4567, Ø 5.6 mm/7 mm TEHA4603, Ø 6.0 mm/3 mm TEHA4605, Ø 6.0 mm/5 mm	 TEHA5683, Ø 6.8 mm/3 mm TEHA5685, Ø 6.8 mm/5 mm							
One-Piece Healing Collars	Eztetic Implants		Trabecular Metal and Tapered Screw-Vent Implants									
	Ø 2.9 mm platform		Ø 3.5 mm platform		Ø 4.5 mm platform		Ø 5.7 mm platform					
Healing Collars 3 mm length shown	 CHCNP31, Ø 3.7 mm/1.5 mm CHCNP33, Ø 3.7 mm/3 mm CHCNP34, Ø 3.7 mm/4.5 mm CHCNP41, Ø 4.5 mm/1.5 mm CHCNP43, Ø 4.5 mm/3 mm CHCNP44, Ø 4.5 mm/4.5 mm	 HC333, Ø 3.5 mm/3 mm HC335, Ø 3.5 mm/5 mm HC343, Ø 4.5 mm/3 mm HC345, Ø 4.5 mm/5 mm HC347, Ø 4.5 mm/7 mm HC353, Ø 5.5 mm/3 mm HC355, Ø 5.5 mm/5 mm	 HC443, Ø 4.5 mm/3 mm HC445, Ø 4.5 mm/5 mm HC453, Ø 5.5 mm/3 mm HC455, Ø 5.5 mm/5 mm HC457, Ø 5.5 mm/7 mm HC463, Ø 6.5 mm/3 mm HC465, Ø 6.5 mm/5 mm	 HC563, Ø 6.5 mm/3 mm HC565, Ø 6.5 mm/5 mm								
Impression Copings	Eztetic Implants		Trabecular Metal and Tapered Screw-Vent Implants									
	Ø 2.9 mm platform		Ø 3.5 mm platform		Ø 4.5 mm platform		Ø 5.7 mm platform					
Indirect Transfer (For closed-tray impressions – includes screw p/n HLTS2)	 CITNP31, Ø 3.7 mm/1.5 mm CITNP33, Ø 3.7 mm/3 mm CITNP34, Ø 3.7 mm/4.5 mm CITNP41, Ø 4.5 mm/1.5 mm CITNP43, Ø 4.5 mm/3 mm CITNP44, Ø 4.5 mm/4.5 mm	 HLT3/3, Ø 3.5 mm HLT3/4, Ø 4.5 mm HLT3/5, Ø 5.5 mm	 HLT4/4, Ø 4.5 mm HLT4/5, Ø 5.5 mm HLT4/6, Ø 6.5 mm	 HLT5/6, Ø 6.5 mm								
Direct Transfer (For open-tray impressions – includes screw p/n DHTS)	 CDTNP31, Ø 3.7 mm/1.5 mm CDTNP33, Ø 3.7 mm/3 mm CDTNP34, Ø 3.7 mm/4.5 mm CDTNP41, Ø 4.5 mm/1.5 mm CDTNP43, Ø 4.5 mm/3 mm CDTNP44, Ø 4.5 mm/4.5 mm	 DHT3/3, Ø 3.5 mm DHT3/4, Ø 4.5 mm DHT3/5, Ø 5.5 mm	 DHT4/4, Ø 4.5 mm DHT4/5, Ø 5.5 mm DHT4/6, Ø 6.5 mm	 HLT5/6 + HLTE, Ø 6.5 mm								
Implant Analogs	 CIANP	 IA3	 IA4	 IA5								
Provisional Restorations	Eztetic Implants		Trabecular Metal and Tapered Screw-Vent Implants									
	Ø 2.9 mm platform		Ø 3.5 mm platform		Ø 4.5 mm platform		Ø 5.7 mm platform					
Temporary Titanium/ Plastic Abutments (Plastic Abutments include retaining screw p/n PTSC and long processing screw p/n DHTS2. for max. 180 days in the case of Plastic Abutments)	 CTANP31, Ø 3.5 mm		 PT341S	 PT344S	 PT341A	 PT344A	 PT451S	 PT454S	 PT451A	 PT454A	 PT561S	 PT564S
Fixture Mount/Transfer ¹ (Incl. replacement retaining screw p/n URS2)			 FMT3	 FMT4	 FMT5							

¹ Supplied with implant.











Zfx™ GenTek™ Ti-Bases and Restorative Components

Crowns and Bridges	Trabecular Metal and Tapered Screw-Vent Implants			
	Eztetic Implants	Ø 3.5 mm platform	Ø 4.5 mm platform	Ø 5.7 mm platform
Ti-Base Engaging Single-Unit, Implant-Level Incl. Retaining and Gluing Screws Recommended Torque 30 Ncm Height 4.7 mm	 ZFX11-ZB-EZ-2947-ES	 ZFX11-ZB-TSV-3547-EL	 ZFX11-ZB-TSV-4547-EL	 ZFX11-ZB-TSV-5747-EL
Ti-Base Non-Engaging Single-Unit, Implant-Level Incl. Retaining and Gluing Screws Recommended Torque 30 Ncm Height 4.7 mm	 ZFX11-ZB-EZ-2947-NES	 ZFX11-ZB-TSV-3547-NEL	 ZFX11-ZB-TSV-4547-NEL	 ZFX11-ZB-TSV-5747-NEL
Retaining Screw For Ti-Base	 ZFX09-ZB-EZ-RS	 ZFX09-ZB-TSV-RS	 ZFX09-ZB-TSV-RS	 ZFX09-ZB-TSV-RS
Ti-Base Non-Engaging Multi-unit, Abutment-Level Incl. Retaining and Gluing Screws Recommended torque 20 Ncm Height 3.9 mm			 ZFX11-ZB-TA-45	
Retaining Screw For Multi-unit Ti-Base			 ZFX09-ZB-TA-RS	
Pre-Abutment Blanks	Trabecular Metal and Tapered Screw-Vent Implants			
	Eztetic Implants	Ø 3.5 mm platform	Ø 4.5 mm platform	Ø 5.7 mm platform
Pre-Abutment Blank (Titanium) Incl. Retaining Screw, 30 Ncm Recommended Torque 30 Ncm	 ZFX18-ZB-EZ-29	 ZFX18-ZB-TSV-35	 ZFX18-ZB-TSV-45	 ZFX18-ZB-TSV-57
Retaining Screw For Ti-Base	 ZFX09-ZB-EZ-RS	 ZFX09-ZB-TSV-RS	 ZFX09-ZB-TSV-RS	 ZFX09-ZB-TSV-RS
Scanbodies	Trabecular Metal and Tapered Screw-Vent Implants			
	Eztetic Implants	Ø 3.5 mm platform	Ø 4.5 mm platform	Ø 5.7 mm platform
Desktop Scanbody, Implant-Level Incl. Screw	 ZFX05D-ZB-EZ-29	 ZFX05D-ZB-TSV-35	 ZFX05D-ZB-TSV-45	 ZFX05D-ZB-TSV-57
Intraoral Scanbodies, Implant-Level Incl. Screws, Double pack Available in heights of 4.0 mm and 7.0 mm	 ZFX05I-ZB-EZ-29	 ZFX05I-ZB-TSV-35	 ZFX05I-ZB-TSV-45	 ZFX05I-ZB-TSV-57
Desktop Scanbody, Abutment Level Incl. Screw			 ZFX05D-ZB-TA-45	
Intraoral Scanbody, Abutment-Level Incl. Screws, Double pack Available in height of 4.0 mm only			 ZFX05I-ZB-TA-45	
Digital Scan Analogs	Trabecular Metal and Tapered Screw-Vent Implants			
	Eztetic Implants	Ø 3.5 mm platform	Ø 4.5 mm platform	Ø 5.7 mm platform
Digital Scan Analog, Implant-Level	 ZFX15-ZB-EZ-29	 ZFX15-ZB-TSV-35	 ZFX15-ZB-TSV-45	 ZFX15-ZB-TSV-57
Placement Tool For Digital Scan Analogs	 ZFX02-M16	 ZFX02-M18	 ZFX02-M18	 ZFX02-M18
Digital Scan Analogs, Abutment Level			 ZFX15-ZB-TA-45	
Placement Tool For Digital Scan Analog Abutment Level			 ZFX02-M18	

Cement-Retained Restorations

Straight	Eztetic Implants		Trabecular Metal and Tapered Screw-Vent Implants		
	Ø 2.9 mm platform	Ø 3.5 mm platform	Ø 4.5 mm platform	Ø 5.7 mm platform	
Hex-Lock® Contour Abutments (Incl. screw p/n MHLAS)	 <p>CANP31S, Ø 3.7 mm / 1.5 mm CANP33S, Ø 3.7 mm / 3 mm CANP34S, Ø 3.7 mm / 4.5 mm CANP41S, Ø 4.5 mm / 1.5 mm CANP43S, Ø 4.5 mm / 3 mm CANP44S, Ø 4.5 mm / 4.5 mm</p>	 <p>ZOA341S, Ø 4.5 mm / 1 mm ZOA342S, Ø 4.5 mm / 2 mm ZOA343S, Ø 4.5 mm / 3 mm</p>	 <p>ZOA441S¹, Ø 4.5 mm / 1 mm ZOA442S², Ø 4.5 mm / 2 mm ZOA443S¹, Ø 4.5 mm / 3 mm ZOA451S, Ø 5.5 mm / 1 mm ZOA452S, Ø 5.5 mm / 2 mm ZOA453S, Ø 5.5 mm / 3 mm</p>	 <p>ZOA561S, Ø 6.5 mm / 1 mm ZOA562S, Ø 6.5 mm / 2 mm ZOA563S, Ø 6.5 mm / 3 mm</p>	
Contour Restorative Kits (Incl. impression cap, analog, provisional and waxing copings. These are also sold individually by extension diameter.)		 <p>ZCRKIT4S Implants Ø 3.7 / 4.1 mm</p>	 <p>ZCRKIT5S Implants Ø 4.7 mm</p>	 <p>ZCRKIT6S Implants Ø 6.0 mm</p>	
Hex-Lock Short Abutments (Incl. screw p/n MHLAS)		 <p>SA341, Ø 4.5 mm / 1 mm SA342, Ø 4.5 mm / 2 mm</p>	 <p>SA451, Ø 5.5 mm / 1 mm SA452, Ø 5.5 mm / 2 mm</p>	 <p>SA561, Ø 6.5 mm / 1 mm SA562, Ø 6.5 mm / 2 mm</p>	
Short Restorative Kits (Incl. impression cap, analog, provisional and waxing copings. These are also sold individually by extension diameter.)		 <p>SARKIT4 Implants Ø 3.7 / 4.1 mm</p>	 <p>SARKIT5 Implants Ø 4.7 mm</p>	 <p>SARKIT6 Implants Ø 6.0 mm</p>	
Hex-Lock Abutments (Incl. screw p/n MHLAS)		 <p>HLA3/3, Ø 3.5 mm HLA3/4, Ø 4.5 mm HLA3/5, Ø 5.5 mm</p>	 <p>HLA4/4, Ø 4.5 mm HLA4/5, Ø 5.5 mm HLA4/6, Ø 6.5 mm</p>	 <p>HLA5/6, Ø 6.5 mm</p>	































¹ Structures for Ø 4.5 mm x Ø 4.5 mm profile are listed but not illustrated. Note: For these items, please use the Contour Restorative Kit ZCRKIT4S (straight) or ZCRKIT4A (angled).

Angled	Eztetic Implants		Trabecular Metal and Tapered Screw-Vent Implants		
	Ø 2.9 mm platform	Ø 3.5 mm platform	Ø 4.5 mm platform	Ø 5.7 mm platform	
Hex-Lock Contour Abutments, 17° (Incl. screw p/n AH20S)	 <p>CANP31A, Ø 3.7 mm / 1.5 mm CANP33A, Ø 3.7 mm / 3 mm CANP41A, Ø 4.5 mm / 1.5 mm CANP43A, Ø 4.5 mm / 3 mm</p>	 <p>ZOA341A, Ø 4.5 mm / 1 mm ZOA342A, Ø 4.5 mm / 2 mm</p>	 <p>ZOA441A¹, Ø 4.5 mm / 1 mm ZOA442A², Ø 4.5 mm / 2 mm ZOA451A, Ø 5.5 mm / 1 mm ZOA452A, Ø 5.5 mm / 2 mm</p>	 <p>ZOA561A, Ø 6.5 mm / 1 mm ZOA562A, Ø 6.5 mm / 2 mm</p>	
Contour Restorative Kits, 17° (Incl. impression cap, analog, provisional and waxing copings. These are also sold individually by extension diameter.)		 <p>ZCRKIT4A Implants Ø 3.7 / 4.1 mm</p>	 <p>ZCRKIT5A Implants Ø 4.7 mm</p>	 <p>ZCRKIT6A Implants Ø 6.0 mm</p>	
Angled Abutment, 20°, for 6 Positions (Incl. screw p/n AH20S)	 <p>C20A3, Ø 3.5 mm</p>	 <p>AH20/4, Ø 4.5 mm</p>	 <p>AH20W/5, Ø 5.5 mm</p>	 <p>A5H20/6, Ø 6.5 mm</p>	




Custom Restorations, Engaging and Non-Engaging





Custom Restorations	Eztetic Implants		Trabecular Metal and Tapered Screw-Vent Implants		
	Ø 2.9 mm platform	Ø 3.5 mm platform	Ø 4.5 mm platform	Ø 5.7 mm platform	
"Cast-To" Gold Abutment, Engaging (Incl. screw p/n MHLAS)		 <p>HLA3G, Ø 3.5 mm</p>	 <p>HLA4G, Ø 4.5 mm</p>	 <p>HLA5G, Ø 5.7 mm</p>	
"Cast-To" Gold Abutment, Non-Engaging (Incl. screw p/n MHLAS)	 <p>CEANP31, Ø 2.9 mm</p>	 <p>NEA3G</p>	 <p>NEA4G</p>		

Replacement Retaining Screws

Screws	Eztetic Implants		Trabecular Metal and Tapered Screw-Vent Implants		
	Ø 2.9 mm platform	Ø 3.5 mm platform	Ø 4.5 mm platform	Ø 5.7 mm platform	
Abutment Screw For Straight Abutments	 <p>CUAS</p>	 <p>MHLAS</p>	 <p>MHLAS</p>	 <p>MHLAS</p>	
Abutment Screw For Angled Abutments	 <p>CUASA Cuff Height, 1.5 mm CUAS Cuff Height, 3 mm</p>	 <p>AH20S</p>	 <p>AH20S</p>	 <p>AH20S</p>	
Replacement Retaining Screw, Direct Transfer	 <p>CASLC</p>	 <p>DHTS</p>	 <p>DHTS</p>	 <p>DHTS</p>	
Replacement Retaining Screw Plastic Temporary Abutments (Long)		 <p>DHTS2</p>	 <p>DHTS2</p>	 <p>DHTS2</p>	
Extension Screw 3.0 mm, Indirect Transfer		 <p>HLTE</p>	 <p>HLTE</p>	 <p>HLTE</p>	
Replacement Retaining Screw, Indirect Transfer, Hex-Lock Abutments (Long)		 <p>HLTS2</p>	 <p>HLTS2</p>	 <p>HLTS2</p>	
Waxing Screw "Cast-To" Gold Abutments, Engaging		 <p>MTWSD</p>	 <p>MTWSD</p>	 <p>MTWSD</p>	
Waxing Screw Tapered Abutment		 <p>SCWS</p>	 <p>SCWS</p>	 <p>SCWS</p>	
Replacement Screw For Copings		 <p>SCTS</p>	 <p>SCTS</p>	 <p>SCTS</p>	

Screw-Retained Restorations




Tapered Abutment	Eztetic Implants		Trabecular Metal and Tapered Screw-Vent Implants		
	Ø 2.9 mm platform	Ø 3.5 mm platform	Ø 4.5 mm platform	Ø 5.7 mm platform	Ø 5.7 mm platform
Tapered Abutments Non-Engaging (Incl. a healing cap p/n TATHC) Note: Abutments do not engage in the internal hex connection. Not for use in single-unit restorations.		 TAC1, 0.75 mm TAC2, 2 mm TAC3, 3 mm TAC4, 4 mm TAC5, 5 mm	 TACW1, 0.75 mm TACW2, 2 mm TACW3, 3 mm TACW4, 4 mm TACW5, 5 mm	 TASC1, 0.75 mm TASC2, 2 mm TASC3, 3 mm TASC4, 4 mm	

Screw-Retained Restorations	Eztetic Implants		Trabecular Metal and Tapered Screw-Vent Implants		
	Ø 2.9 mm platform	Ø 3.5 mm platform	Ø 4.5 mm platform	Ø 5.7 mm platform	Ø 5.7 mm platform
Angled Tapered Abutments 15° Engaging in implant, not for use in single-unit restorations. (Incl. preassembled retaining screw)		 15AT312, 1 mm / 2 mm 15AT323, 2 mm / 3 mm 15AT334, 3 mm / 4 mm	 15AT412, 1 mm / 2 mm 15AT423, 2 mm / 3 mm 15AT434, 3 mm / 4 mm		
Angled Tapered Abutments 30° Engaging in implant, not for use in single-unit restorations. (Incl. preassembled retaining screw)		 30AT324, 2 mm / 4 mm 30AT335, 3 mm / 5 mm 30AT346, 4 mm / 6 mm	 30AT424, 2 mm / 4 mm 30AT435, 3 mm / 5 mm 30AT446, 4 mm / 6 mm		

Tapered Abutment Components

								
ACTIT	ACTDT	ACTR	ACTGC	ACTP	TTC5	ACTT	TGC3	TGC5
Tapered Abutment Indirect and Direct Transfers	Tapered Abutment Replica (Represents Tapered Abutment Attached to Implant)	Tapered Abutment Gold Coping with Screw (p/n SCTS)	Plastic Castable Coping with Screw (p/n SCTS)	Tapered Abutment Titanium Bar Coping with Screw (p/n SCTS)	Titanium Temporary Coping with Screw (p/n SCTS)	Tapered Abutment Bar Gold Copings with Screw		








Overdenture Restorations





LOCATOR R-Tx® Abutments	Eztetic Implants		Trabecular Metal and Tapered Screw-Vent Implants		
	Ø 2.9 mm platform	Ø 3.5 mm platform	Ø 4.5 mm platform	Ø 5.7 mm platform	Ø 5.7 mm platform
LOCATOR R-Tx Abutments		 30200-00, 0.5 mm 30200-01, 1 mm 30200-02, 2 mm 30200-03, 3 mm 30200-04, 4 mm 30200-05, 5 mm 30200-06, 6 mm	 30201-00, 0.5 mm 30201-01, 1 mm 30201-02, 2 mm 30201-03, 3 mm 30201-04, 4 mm 30201-05, 5 mm 30201-06, 6 mm	 30202-01, 1 mm 30202-02, 2 mm 30202-03, 3 mm 30202-04, 4 mm 30202-05, 5 mm 30202-06, 6 mm	

LOCATOR R-Tx Lab Processing Components









					
30017-01	30014-01	30015-01	30016-01	30018-01	30013-01
Impression Coping	0.35 mmD Abutment Analog	4.0 mmD Abutment Analog	5.0 mmD Abutment Analog	Processing Spacer	Denture Attachment Processing Assembly

LOCATOR R-Tx Retention Inserts





						
08514	30001-01	30002-01	30003-01	30004-01	30012-01	30021-01
Block-out Spacer	Zero Retention Insert	Low Retention Insert	Medium Retention Insert	High Retention Insert	Processing Insert	Retention Insert Tool

LOCATOR® Abutments	Eztetic Implants		Trabecular Metal and Tapered Screw-Vent Implants		
	Ø 2.9 mm platform	Ø 3.5 mm platform	Ø 4.5 mm platform	Ø 5.7 mm platform	Ø 5.7 mm platform
LOCATOR Abutments	 CLOCNP-1, 1 mm CLOCNP-2, 2 mm CLOCNP-3, 3 mm CLOCNP-4, 4 mm CLOCNP-5, 5 mm CLOCNP-6, 6 mm	 TLOC3/0, 0.5 mm TLOC3/1, 1 mm TLOC3/2, 2.5 mm TLOC3/3, 3.5 mm TLOC3/4, 4.5 mm TLOC3/5, 5.5 mm TLOC3/6, 6.5 mm	 TLOC4/0, 0.5 mm TLOC4/1, 1 mm TLOC4/2, 2 mm TLOC4/3, 3 mm TLOC4/4, 4 mm TLOC4/5, 5 mm TLOC4/6, 6 mm	 TLOC5/1, 1 mm TLOC5/2, 2 mm TLOC5/3, 3 mm TLOC5/4, 4 mm	

LOCATOR Attachment Components

							
LOCIMP	LOCAN/4	LOCP	LOCBF2	LOCTB	LOCLBF	LOCMP2	LOCMPX2
Transfer Cap	Abutment Analogue	Parallel Coping	Bar Female	"Cast-To" Bar Female	Laser Bar Female	Male Processing Package	Male Processing Package Extended Range




Note: Replacement Males are available. See product Item.

Ball Abutments	Eztetic Implants		Trabecular Metal and Tapered Screw-Vent Implants		
	Ø 2.9 mm platform	Ø 3.5 mm platform	Ø 4.5 mm platform	Ø 5.7 mm platform	Ø 5.7 mm platform
Ball Abutments with Metal Housing and Nylon Liners	 CBANP2, Ø 2.9 mm / 2 mm CBANP4, Ø 2.9 mm / 4 mm CBANP6, Ø 2.9 mm / 6 mm	 BAC2, 2 mm BAC4, 4 mm BAC6, 6 mm	 BACW2, 2 mm BACW4, 4 mm BACW6, 6 mm	 BA5C2, 2 mm BA5C4, 4 mm	

Transfer Components for Ball Abutments

					
BAT	CAT	BAR	CAH	CAN	CAN-G
Ball Abutment Transfer Cap	Transfer Cap	Ball Abutment Analogue	Metal Females	Nylon Liner	Nylon Liner for Fixed Retention

Prosthetic Tools and Torque Wrenches

					
LOCCT2	30021-01	LOCTW	ABTH		
LOCATOR Tool, 3-piece	Tool for R-Tx LOCATOR	LOCATOR Torque Wrench Insert	Abutment Holder		
					
TWR	TW1.25	TW1.25L	HXGR1.25	HX1.25D	TLRT2
Adjustable Torque Wrench 10-35 Ncm	Torque Wrench Insert Hex Driver, Short, 1.25 mmD, 17 mm Long	Torque Wrench Insert Hex Driver, Long, 1.25 mmD, for TWR, 22 mm Long	Hex Driver, Short 22 mm Long	Latch-Lock Hex Driver 23 mm long	Removal Tool for Internal Hex Implant Abutments
					
			HXLGR1.25	HXL1.25D	CLRT2
			Hex Driver, Long 30 mm Long	Latch-Lock Hex Driver 30 mm long	Removal Tool for Eztetic Abutments



Contact us at 1-800-342-5454 or visit
zimmerbiometdental.com

Zimmer Biomet Dental
Global Headquarters
4555 Riverside Drive
Palm Beach Gardens, FL 33410
Tel: +1-561-776-6700
Fax: +1-561-776-1272

Unless otherwise indicated, as referenced herein, all trademarks are the property of Zimmer Biomet; and all products are manufactured by one or more of the dental subsidiaries of Zimmer Biomet Holdings, Inc. and marketed and distributed by Zimmer Biomet Dental and its authorized marketing partners. LOCATOR and LOCATOR R-Tx are trademarks of Zest IP Holdings, LLC. LOCATOR Abutment and LOCATOR R-Tx Attachment Systems are manufactured by Zest IP Holdings, LLC. For additional product information, please refer to the individual product labeling or instructions for use. Product clearance and availability may be limited to certain countries/regions. This material is intended for clinicians only and does not comprise medical advice or recommendations. Distribution to any other recipient is prohibited. This material may not be copied or reprinted without the express written consent of Zimmer Biomet Dental. ZB1194 REV A 07/21 ©2021 Zimmer Biomet. All rights reserved.

